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Parameter	Function
<b>--outSAMattributes spatial</b>	Set to turn on spatial BAM file format mode.
<b>--outSAMtype BAM SortedByCoordinate</b>	(STAR option) Set output BAM file sorted by coordinate.
<b>--genomeDir</b>	(STAR option) Path to the directory where the genome indices are stored.
<b>--runThreadN</b>	(STAR option; defaults to 1) Set the number of threads to be used.
<b>--outFileNamePrefix</b>	(STAR option) Custom output file prefix.
<b>--sysShell /bin/bash</b>	(STAR option) Path to the shell binary.
<b>--stParaFile</b>	(Required) Name of a parameters file defines CID mapping options. Options are specified in Table 2-2.
<b>--readNameSeparator \ " \ "</b>	(STAR option) Character(s) separating the part of the read names that will be trimmed in output.
<b>--limitBAMsortRAM</b>	(STAR option) Maximum available RAM (bytes) for sorting BAM.
<b>--limitOutSJcollapsed</b>	(STAR option) Max number of collapsed junctions.
<b>--limitIObufferSize</b>	(STAR option) Max available buffers size (bytes) for input/output, per thread.
<b>--outBAMsortingBinsN</b>	(STAR option; defaults to 50) number of genome bins for coordinate-sorting. If the read2 FASTQ file size is greater than 200, it is better to set this value to 100.

表 2-2 mapping --stParaFile 命令行参数

Parameter	Function
<b>in</b>	(Required) Path to the STOmics Chip Mask file.
<b>in1</b>	(Required) Path to the FASTQ file. If PE sequencing, then specify the path to the FASTQ file of read1 here.
<b>in2</b>	(Optional) Path to the FASTQ file of read2. Only valid for PE sequencing.
<b>encodeRule</b>	(Required) Encoding rule for the four bases. ACTG stands for A->0, C->1, T->2, G->3.
<b>out</b>	(Optional) Set output file prefix.
<b>action</b>	(Required; defaults to 1) Action number. Valid options: 1=CID stat, 2=CID overlap, 3=get CID position map, 4=map CID to slide, or 5=merge CID list.
<b>barcodeReadsCount</b>	(Required) Mapped CID list file with reads counts for each CID.
<b>platform</b>	(Optional) Sequencing platform. Valid options: SEQ500, T1, or T10.
<b>barcodeStart</b>	(Required; defaults to 0) CID start position.
<b>barcodeLen</b>	(Required; defaults to 25) CID length.
<b>umiStart</b>	(Required; defaults to 25) MID start position.
<b>umiLen</b>	(Required; defaults to 10) MID length.
<b>umiRead</b>	(Required; defaults to 1) Declare the read contains MID.
<b>mismatch</b>	(Required; defaults to 0) Max mismatch tolerant.
<b>bcNum</b>	(Optional) CID count. Please check 2.9.1 CIDcount for more information.

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Parameter	Function
<b>-i</b>	(Required) Mapped CID list files with reads counts for each CID.
<b>--out</b>	(Required) Mapped CID list file which merges all input files.
<b>--action</b>	(Required; defaults to 1) Action number. Valid options: 1=map CID to slide, 2=merge CID list, 3=mask format change, or 4=mask merge.

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```

$ mkdir /path/to/multi_lane_output/01.merge
$ singularity exec SAW_v4.1.0.sif merge \
    -i /path/to/multi_lane_output/00.mapping/{lane1}.barcodeReadsCount.
[*] txt,/path/to/multi_lane_output/00.mapping/{lane2}.barcodeReadsCount.txt \
    --out /path/to/multi_lane_output/01.merge/{SN}.barcodeReadsCount.txt
    --action 2

```

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Parameter	Function
<b>-i</b>	(Required) mapping output BAM file. Separate multiple files by comma.
<b>-o</b>	(Required) Set the count output BAM file name.
<b>-a</b>	(Required) Gene annotation file.
<b>-s</b>	(Required) Set the count output statistical summary report file name.
<b>-e</b>	(Required) Set the count output gene expression file name.
<b>--umi_len</b>	(Required; defaults to 10) MID length.
<b>--sn</b>	(Required) STOmics Chip serial number (SN).
<b>-c</b>	(Optional) CPU core number to use.
<b>--save_lq</b>	(Optional) Save low quality reads if set.
<b>--save_dup</b>	(Optional) Save duplicate reads if set.
<b>--umi_on</b>	(Optional) Correct MID if set.
<b>--sat_file</b>	(Optional) Set the saturation sampling file name which is prepared for sequencing saturation (requires <b>--umi_on</b> ).
<b>-m</b>	(Optional; defaults to detected) Set available memory (GB).

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```

$ mkdir -p /path/to/output/02.count
$ geneExp=/path/to/output/02.count/{SN}.raw.gef
$ saturationSamplingFile=/path/to/output/02.count/raw_barcode_gene_exp.txt
$ singularity exec SAW_v4.1.0.sif count \
    -i /path/to/output/00.mapping/{lane}.Aligned.sortedByCoord.out.bam \
    -o /path/to/output/02.count/{SN}.Aligned.sortedByCoord.out.merge.q10.dedup.
[*] target.bam \
    -a /path/to/genomeDir/genes.gtf \
    -s /path/to/output/02.count/{SN}.Aligned.sortedByCoord.out.merge.q10.dedup.
[*] target.bam.summary.stat \
    -e ${geneExp} \
    --umi_len 10 \
    --sat_file ${saturationSamplingFile} \
    --sn {SN} \
    --umi_on \
    --save_lq \
    --save_dup \
    -c 8 \
    -m 128

```

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Parameter	Function
<b>-i</b>	(Required) ImageQC processed staining image TAR.GZ file.
<b>-c</b>	(Required) ImageQC JSON file.
<b>-v</b>	(Required) count output gene expression matrix GEF file.
<b>-o</b>	(Required) Path to the directory where to store the register outputs.

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```

$ image=/path/to/data/image
$ image4register=$(find ${image} -maxdepth 1 -name {SN}*.tar.gz | head -1)
$ imageQC=$(find ${image} -maxdepth 1 -name {SN}*.json | head -1)
$ mkdir -p /path/to/output/03.register
$ singularity exec SAW_v4.1.0.sif register \
    -i ${image4register} \
    -c ${imageQC} \
    -v /path/to/output/02.count/{SN}.raw.gef \
    -o /path/to/output/03.register

```

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Parameter	Function
<b>--dnbfile</b>	(Required) Mapped CID list file with reads counts for each CID. Use the merged mapped CID with reads count file if multiple pairs of FASTQ files are involved in the analysis.
<b>-i</b>	(Required) count output gene expression matrix GEF file.
<b>-o</b>	(Required) Path to the directory where to store the <b>tissueCut</b> outputs.
<b>-s</b>	(Optional) Path to the directory where the <b>register</b> outputs are stored. Only valid when <b>register</b> has performed.
<b>-t</b>	(Required) Run for tissue segmentation. Valid options: <b>tissue</b> .
<b>--snId</b>	(Required) STOmics Chip serial number (SN).
<b>--platform</b>	(Required, default to T1) Sequencing platform.

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```
$ mkdir -p /path/to/output/04.tissuecut
$ singularity exec SAW_v4.1.0.sif tissueCut \
  --dnbfile /path/to/output/01.merge/{SN}.barcodeReadsCount.txt \
  -i /path/to/output/02.count/{SN}.raw.gef \
  -o /path/to/output/04.tissuecut \
  -s /path/to/output/03.register/7_result \
  -t tissue \
  --snId {SN} \
  --platform T10
```

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```
$ mkdir -p /path/to/output/04.tissuecut
$ singularity exec SAW_v4.1.0.sif tissueCut \
  --dnbfile /path/to/output/01.merge/{SN}.barcodeReadsCount.txt \
  -i /path/to/output/02.count/{SN}.raw.gef \
  -o /path/to/output/04.tissuecut \
  -t tissue \
  --snId {SN} \
  --platform T10
```

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表 2-7 spatialCluster 命令行参数

Parameter	Function
<b>-i</b>	(Required) tissuCut output GEF file for the tissue coverage area.
<b>-o</b>	(Required) Output path for the clustering result in H5AD format.
<b>-s</b>	(Required) Bin size, recommend to use 200.

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```
$ mkdir -p /path/to/output/05.spatialcluster
$ singularity exec SAW_v4.1.0.sif spatialCluster \
  -i /path/to/output/04.tissuecut/{SN}.tissue.gef \
  -o /path/to/output/05.spatialcluster/{SN}.spatial.cluster.h5ad \
  -s 200
```

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```
$ ll -Rht /path/to/output/05.spatialcluster
-rw-rw-r-- 2 ubuntu ubuntu 231M Apr 13 22:40 SN.spatial.cluster.h5ad
```

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Parameter	Function
<b>-i</b>	(Required) count output saturation sampling file.
<b>--tissue</b>	(Required) tissuCut output GEF file for the tissue coverage area.
<b>-o</b>	(Required) Path to the directory where to store the saturation outputs.
<b>--bcstat</b>	(Required) mapping output statistical report of CID mapping. Separate multiple files by comma.
<b>--summary</b>	(Required) count output statistical report of annotation.

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```
$ mkdir -p /path/to/output/06.saturation
$ singularity exec SAW_v4.1.0.sif saturation \
  -i /path/to/output/02.count/raw_barcode_gene_exp.txt \
  --tissue /path/to/output/04.tissuecut/{SN}.tissue.gef \
  -o /path/to/output/06.saturation \
  --bcstat /path/to/output/00.mapping/{lane}_barcodeMap.stat \
  --summary /path/to/output/02.count/{SN}.Aligned.sortedByCoord.out.merge.q10.
[*] dedup.target.bam.summary.stat
```

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Parameter	Function
<b>-m</b>	(Required) Statistical report of CID mapping. Separate multiple files by comma.
<b>-a</b>	(Required) Statistical report of STAR alignment. Separate multiple files by comma.
<b>-g</b>	(Required) Statistical report of annotation.
<b>-l</b>	(Required) Statistical report of tissue-covered region.
<b>-n</b>	(Required) tissueCut output GEF file.
<b>-b</b>	(Required) tissueCut output bin 200 scatter plot.
<b>-v</b>	(Required) tissueCut output bin 200 violin plots.
<b>-i</b>	(Optional) The image pyramid RPI file.
<b>-d</b>	(Required) spatialCluster output H5AD file.
<b>-t</b>	(Required) saturation output bin 200 sequence saturation plot.
<b>-r standard_version</b>	(Required) Set to specifying report version.
<b>-s</b>	(Required) The STOmics Chip serial number.
<b>--pipelineVersion</b>	(Required) Set to specifying analysis pipeline version.
<b>-o</b>	(Required) The directory to store outputs.

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```

$ mkdir -p /path/to/output/07.report
$ singularity exec SAW_v4.1.0.sif report \
  -m /path/to/output/00.mapping/{lane}_barcodeMap.stat \
  -a /path/to/output/00.mapping/{lane}.Log.final.out \
  -g /path/to/output/02.count/{SN}.Aligned.sortedByCoord.out.merge.q10.dedup.
[*]target.bam.summary.stat \
  -l /path/to/output/04.tissuecut/tissuecut.stat \
  -n /path/to/output/04.tissuecut/{SN}.gef \
  -d /path/to/output/05.spatialcluster/{SN}.spatial.cluster.h5ad \
  -t /path/to/output/06.saturation/plot_200x200_saturation.png \
  -b /path/to/output/04.tissuecut/tissue_fig/scatter_200x200_MID_gene_counts.png \
  -v /path/to/output/04.tissuecut/tissue_fig/violin_200x200_MID_gene.png \
  -r standard_version \
  -i /path/to/output/04.tissuecut/tissue_fig/{SN}.ssDNA.rpi \
  -s {SN} \
  --pipelineVersion SAW_v4.1.0 \
  -o /path/to/output/07.report

```

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```
$ singularity exec SAW_v4.1.0.sif CIDcount \
  -i /path/to/data/{SN}.barcodeToPos.h5 \ ## STOmics Chip Mask file path
  -s {speciesName} \ ## a string of species name
  -g {genomeSize} ## Genome file size in GB, can be acquired by "ls -l --block-
size=GB ${Genome file of the species after STAR indexing}"
```

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```
...  
$ singularity exec SAW_v4.1.0.sif gefTools bgef \  ## complete GEF that only contains  
bin1 geneExp group to a whole GEF, you may specify the bin size you need using "-b".  
Separate multiple bin size with comma  
-i /path/to/output/02.count/{SN}.tissue.gef \  
-o {SN}.tissue.complete.gef \  
-b 1,20,50,100
```

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```
...  
$ singularity exec SAW_v4.1.0.sif gefTools bgef \  ## convert GEM to GEF in specific  
bin size. Separate multiple bin sizes with comma  
-i {SN}.gem \  
-o {SN}.gef \  
-b 1,20,50
```

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```
$ cat /path/to/output/00.mapping/E100026571_L01_trim_read_barcodeMap.stat
...
getBarcodePositionMap_uniqBarcodeTypes: 645784920
total_reads: 1002214171
fixed_sequence_containing_reads: 0 0.00%
pass_filter_reads: 1002214171
mapped_reads: 845170516 84.33%
reads_with_adapter: 8137401 0.81%
reads_with_dnb: 45284608 4.52%
barcode_exactly0overlap_reads: 698287595 69.67%
barcode_mis0overlap_reads: 146882921 14.66%
barcode_withN_reads: 0 0.00%
Q10_bases_in_barcode: 99.54%
Q20_bases_in_barcode: 97.49%
Q30_bases_in_barcode: 91.74%
Q10_bases_in_umi: 99.26%
Q20_bases_in_umi: 96.32%
Q30_bases_in_umi: 89.45%
Q10_bases_in_seq: 99.47%
Q20_bases_in_seq: 97.12%
Q30_bases_in_seq: 91.08%
umi_filter_reads: 8451821 0.84%
umi_with_N_reads: 13355 0.00%
umi_with_polyA_reads: 13044 0.00%
umi_with_low_quality_base_reads: 8425422 0.84%
mapped_dnbs: 78023582
...
```

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```

$ cat /path/to/output/00.mapping/E100026571_L01_trim_read.Log.final.out
...
          Number of input reads | 783296686
    Average input read length | 100
          UNIQUE READS:
    Uniquely mapped reads number | 624859097
      Uniquely mapped reads % | 79.77%
        Average mapped length | 98.02
      Number of splices: Total | 68091842
Number of splices: Annotated (sjdb) | 65321302
      Number of splices: GT/AG | 66183758
      Number of splices: GC/AG | 462810
      Number of splices: AT/AC | 43237
Number of splices: Non-canonical | 1402037
      Mismatch rate per base, % | 0.58%
        Deletion rate per base | 0.07%
        Deletion average length | 3.87
        Insertion rate per base | 0.04%
        Insertion average length | 1.26
          MULTI-MAPPING READS:
    Number of reads mapped to multiple loci | 81769070
      % of reads mapped to multiple loci | 10.44%
    Number of reads mapped to too many loci | 2071180
      % of reads mapped to too many loci | 0.26%
          UNMAPPED READS:
    Number of reads unmapped: too many mismatches | 0
      % of reads unmapped: too many mismatches | 0.00%
    Number of reads unmapped: too short | 73977028
      % of reads unmapped: too short | 9.44%
    Number of reads unmapped: other | 620311
      % of reads unmapped: other | 0.08%
          CHIMERIC READS:
      Number of chimeric reads | 0
        % of chimeric reads | 0.00%

```

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```
$ cat /path/to/output/04.tissuecut/tissuecut.stat
##### Tissue Statistic Analysis with Stain Image #####
Contour_area: 88336939
Number_of_DNB_under_tissue: 35918358
Ratio: 40.66%
Total_gene_type: 24296
MID_counts: 87120944
Fraction_MID_in_spots_under_tissue: 83.20%
Reads_under_tissue: 657982650
Fraction_reads_in_spots_under_tissue: 77.85%

binSize=1
Mean_reads_per_spot: 18.319
Median_reads_per_spot: 12.0
Mean_gene_type_per_spot: 1.7
Median_gene_type_per_spot: 1.0
Mean_Umi_per_spot: 2.426
Median_Umi_per_spot: 2.0

binSize=50
Mean_reads_per_spot: 18371.193
Median_reads_per_spot: 16481.5
Mean_gene_type_per_spot: 1131.942
Median_gene_type_per_spot: 1098.0
Mean_Umi_per_spot: 2432.459
Median_Umi_per_spot: 2246.0

binSize=100
Mean_reads_per_spot: 72369.407
Median_reads_per_spot: 65597.0
Mean_gene_type_per_spot: 3043.033
Median_gene_type_per_spot: 3038.0
Mean_Umi_per_spot: 9582.154
Median_Umi_per_spot: 8820.5

binSize=150
Mean_reads_per_spot: 160483.573
Median_reads_per_spot: 146299.5
Mean_gene_type_per_spot: 4843.486
Median_gene_type_per_spot: 4969.5
Mean_Umi_per_spot: 21249.011
Median_Umi_per_spot: 19677.0

binSize=200
Mean_reads_per_spot: 281189.167
Median_reads_per_spot: 258306.0
Mean_gene_type_per_spot: 6342.139
Median_gene_type_per_spot: 6640.5
Mean_Umi_per_spot: 37231.173
Median_Umi_per_spot: 34612.5
```

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


```
$ head /path/to/output/new_final_result.json
{
  "version": "version_v2",
  "1.Filter_and_Map": {
    "1.1.Adapter_Filter": [
      {
        "Sample_id": "E100026571_L01_trim_read",
        "getCIDPositionMap_uniqCIDTypes": "645.78M",
        "total_reads": "1.0G",
        "fixed_sequence_contianing_reads": "0.0(0.00%)",
        "mapped_reads": "845.17M(84.33%)",
```

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1. BGIResearch/SAW. Accessed October 13, 2021. <https://github.com/BGIResearch/SAW>
2. Chen A, Liao S, Cheng M, et al. Large field of view-spatially resolved transcriptomics at nanoscale resolution Short title: DNA nanoball stereo-sequencing. *bioRxiv*. Published online January 24, 2021:2021.01.17.427004. doi:10.1101/2021.01.17.427004
3. Cock PJA, Fields CJ, Goto N, Heuer ML, Rice PM. The Sanger FASTQ file format for sequences with quality scores, and the Solexa/Illumina FASTQ variants. *Nucleic Acids Res*. 2009;38(6):1767-1771. doi:10.1093/nar/gkp1137
4. Archives SR, Sra T, Nucleotide I, et al. File Format Guide 1. Published online 2009:1-11. Accessed May 21, 2021. <https://www.ncbi.nlm.nih.gov/sra/docs/submitformats/>
5. Merkel D. Docker: lightweight Linux containers for consistent development and deployment. *Linux J*. 2014;2014(239):2. Accessed October 15, 2021. <https://www.linuxjournal.com/content/docker-lightweight-linux-containers-consistent-development-and-deployment>
6. Kurtzer GM, Sochat V, Bauer MW. Singularity: Scientific containers for mobility of compute. *PLoS One*. 2017;12(5):e0177459. doi:10.1371/journal.pone.0177459
7. *Sequence Alignment/Map Format Specification*.; 2021. Accessed May 21, 2021. <https://github.com/samtools/hts-specs>.
8. Dobin A, Davis CA, Schlesinger F, et al. STAR: Ultrafast universal RNA-seq aligner. *Bioinformatics*. 2013;29(1):15-21. doi:10.1093/bioinformatics/bts635
9. Ensembl. GFF/GTF File Format. Published 2020. Accessed May 27, 2021. <http://www.ensembl.org/info/website/upload/gff.html?redirect=no>
10. GFF2 - GMOD. Accessed May 27, 2021. <http://gmod.org/wiki/GFF2>
11. GitHub - BGIResearch/stereopy: A toolkit of spatial transcriptomic analysis. Accessed July 4, 2021. <https://github.com/BGIResearch/stereopy>
12. BGIResearch/geftools: Tools for manipulating GEFs. Accessed April 7, 2022. <https://github.com/BGIResearch/geftools>
13. BGIResearch/gefpy: gef io, draw out from stereopy. Accessed April 7, 2022. <https://github.com/BGIResearch/gefpy>

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©2022 深圳华大生命科学研究院保留所有权利。